


A Simple Field Assay for Insect-Borne Diseases

Alex Flores,¹ LuAnn Scott,^{1,2} Sanchita Bhadra,³ Andrew D. Ellington,³ Holly A. Wichman^{1,2}

¹Department of Biological Sciences, University of Idaho, Moscow, ID 83844 ²Center for Modeling Complex Interactions, University of Idaho, Moscow, ID 83844 ³University of Texas at Austin, Austin, TX 78712





Background

- Insect-borne diseases continue to expand across a variety of geographical locations.
- Expensive instrumentation for PCR makes monitoring the spread of these diseases difficult.





Figure 1: Global Spread of Zika Virus, 2007-2015. Vectors: Aedes aegypti & A. albopictus

Although Aedes aegypti is the primary vector of many disease-causing viruses (malaria, yellow fever, dengue, chikungunya, and zika), A. albopictus is providing a vector for new viral strains and has spread to all continents in the last 50 years.








Figure 3: Global Spread of Dengue Virus 1, 1949-2013 Vector: Aedes aegypti & A. albopictus

Methods: *Drosophila* Model System Results

- Large populations
- Easy to grow
- Infected by injection (positive controls) or 'naturally' by oral ingestion (test flies).
- 2 RNA viruses grown in cell culture are used: Drosophila C Virus (DCV) Drosophila X Virus (DXV)
- 10 flies sampled each day over the course of infection and assayed individually.


Methods: LAMP-OSD Assay

- Carefully design primers and probes.
- Amplify viral nucleic acids (loop-mediated isothermal amplification, LAMP).
- Transduce amplicons to visible fluorescence via unique toehold-mediated strand displacement reporters (OSD probes).

F1cF2cF3c B3cB2cB1c F1 F2 F3 Hammerhead initiate


Visualize with blue light and record readout with cell phone camera.

- Both viruses were successfully amplified & visualized using ~10² copies of stock.
- Orally infected flies yield various densities of fluorescence.



Raw Data – DCV day 2

- All dead flies were assayed plus enough live flies to total 10 per
- **Negative control without primers** was run for each sampled fly.
- Fluorescent density of mock fly was used to determine cut-off for negative/positive result.



Viral Dynamics Summary

- DCV flies sampled over 10 days
 - Proportion of infected flies decreased after 6 days postinfection.
- DXV flies sampled over 30 days
- Proportion of infected flies decreased after 20 days postinfection.

Conclusions

- LAMP-OSD assays can identify the presence of virus in individual flies.
- Field studies of virus prevalence in arthropod populations can be carried out with minimal equipment: an incubator, a blue diode, and a cell phone.
- Course of infection data for both DCV and DXV is consistent with results from RT-PCR.

In Progress

- **Experiments comparing LAMP-OSD and RT-**PCR in detection of virus.
- Analysis of changes in viral dynamics in the presence of the *Drosophila* bacterial endosymbiont Wolbachia.



- 1. Develop a genetic method that can type hundreds of insects quickly under field conditions.
- 2. Test the method in a *Drosophila* system to monitor infection and transmission in a controlled setting.
 - Develop smartphone-imaged LAMP-OSD assays to detect 2 Drosophila RNA viruses (DCV & DXV).
 - Determine limits of detection of fly viruses.
 - Test assay's capability to monitor a lab-based course of infection.

Acknowledgements

This work is made possible by a BEACON grant from the National Science Foundation. DCV and DXV stock viruses were provided by Tanya Miura. Drosophila stocks are maintained by the Christine Parent lab.

Advice on fly husbandry and infection was provided by Ashley De Aguero & Andrea Gonzalez.

Insect images from Petey Mosquitey and MegaCatch.com.